

Need a new gel tank?

Get one for free with only five boxes of precast gels



Simply quote promo code **R5GEL3** when ordering any five of the qualifying gels to receive a free Mini Gel Tank. Offer available until June 30th 2019*

Novex WedgeWell Tris-Glycine gels	Bolt gels	NuPAGE SDS-PAGE gels
<p>Discover a whole new depth to Tris-Glycine gels with 60 μL load volume</p> <p>Invitrogen™ Novex™ WedgeWell™ Tris-Glycine gels are precast mini gels for separation of proteins by polyacrylamide gel electrophoresis allowing the use of Laemmli sample and running buffers. These gels enable high quality performance and separation of a wide range of proteins into well-resolved bands.</p> <ul style="list-style-type: none"> • Wedge-shaped wells — easily load up to 60 μL of sample per well • Improved shelf life — store gels for up to 12 months at 4°C • Fast run conditions — quickly separate proteins using constant voltage in less than 60 minutes <p>Find out more at thermofisher.com/novexwedge</p>	<p>Designed to deliver unmatched western blot results at an affordable price</p> <p>Precast polyacrylamide gels designed for optimal separation of a broad molecular weight range of proteins under denaturing conditions during gel electrophoresis. Invitrogen™ Bolt™ gels are ideal for western blot transfer and analysis, and any other techniques where protein integrity is crucial.</p> <ul style="list-style-type: none"> • Superior band quality and band volume — Invitrogen™ Bolt™ Bis-Tris Plus chemistry delivers sharp, straight, uniform bands with higher band volume • Better protein resolution due to 10% greater resolving distance and optimised gradient format • Preserved protein integrity with neutral-pH formulation that minimises protein modifications <p>Find out more at thermofisher.com/bolt</p>	<p>The gold standard for superior protein separation</p> <p>Revolutionary high-performance polyacrylamide gel electrophoresis system that simulates the denaturing conditions of the traditional Laemmli system (Invitrogen™ Tris-Glycine SDS-PAGE gels). Invitrogen™ NuPAGE™ SDS-PAGE gels use a unique buffer formulation that maintains a low operating pH during electrophoresis. This helps eliminate the “smiles” and poor resolution seen with Tris-Glycine SDS-PAGE gels.</p> <ul style="list-style-type: none"> • Superior protein band resolution and stability • Faster sample run times—35–50 minutes • Longer product shelf life—16 months from the date of manufacture <p>Find out more at thermofisher.com/nupage</p>

Qualifying Novex WedgeWell Tris-Glycine mini gels

Description	Cat. no.
Novex WedgeWell 6% Tris-Glycine Mini Gels, 10-well	XP00060BOX
Novex WedgeWell 6% Tris-Glycine Mini Gels, 12-well	XP00062BOX
Novex WedgeWell 6% Tris-Glycine Mini Gels, 15-well	XP00065BOX
Novex WedgeWell 8% Tris-Glycine Mini Gels, 10-well	XP00080BOX
Novex WedgeWell 8% Tris-Glycine Mini Gels, 12-well	XP00082BOX
Novex WedgeWell 8% Tris-Glycine Mini Gels, 15-well	XP00085BOX
Novex WedgeWell 10% Tris-Glycine Mini Gels, 10-well	XP00100BOX
Novex WedgeWell 10% Tris-Glycine Mini Gels, 12-well	XP00102BOX
Novex WedgeWell 10% Tris-Glycine Mini Gels, 15-well	XP00105BOX
Novex WedgeWell 12% Tris-Glycine Mini Gels, 10-well	XP00120BOX
Novex WedgeWell 12% Tris-Glycine Mini Gels, 12-well	XP00122BOX
Novex WedgeWell 12% Tris-Glycine Mini Gels, 15-well	XP00125BOX
Novex WedgeWell 14% Tris-Glycine Mini Gels, 10-well	XP00140BOX
Novex WedgeWell 14% Tris-Glycine Mini Gels, 12-well	XP00142BOX
Novex WedgeWell 14% Tris-Glycine Mini Gels, 15-well	XP00145BOX
Novex WedgeWell 16% Tris-Glycine Mini Gels, 10-well	XP00160BOX
Novex WedgeWell 16% Tris-Glycine Mini Gels, 12-well	XP00162BOX
Novex WedgeWell 16% Tris-Glycine Mini Gels, 15-well	XP00165BOX
Novex WedgeWell 4-12% Tris-Glycine Mini Gels, 10-well	XP04120BOX
Novex WedgeWell 4-12% Tris-Glycine Mini Gels, 12-well	XP04122BOX
Novex WedgeWell 4-12% Tris-Glycine Mini Gels, 15-well	XP04125BOX
Novex WedgeWell 4-20% Tris-Glycine Mini Gels, 10-well	XP04200BOX

Description	Cat. no.
Novex WedgeWell 4-20% Tris-Glycine Mini Gels, 12-well	XP04202BOX
Novex WedgeWell 4-20% Tris-Glycine Mini Gels, 15-well	XP04205BOX
Novex WedgeWell 8-16% Tris-Glycine Mini Gels, 10-well	XP08160BOX
Novex WedgeWell 8-16% Tris-Glycine Mini Gels, 12-well	XP08162BOX
Novex WedgeWell 8-16% Tris-Glycine Mini Gels, 15-well	XP08165BOX
Novex WedgeWell 10-20% Tris-Glycine Mini Gels, 10-well	XP10200BOX
Novex WedgeWell 10-20% Tris-Glycine Mini Gels, 12-well	XP10202BOX
Novex WedgeWell 10-20% Tris-Glycine Mini Gels, 15-well	XP10205BOX

Qualifying Bolt Bis-Tris Plus mini gels

Description	Cat. no.
Bolt 8% Bis-Tris Plus Gel, 10 Well	NW00080BOX
Bolt 8% Bis-Tris Plus Gel, 12 Well	NW00082BOX
Bolt 8% Bis-Tris Plus Gel, 15 Well	NW00085BOX
Bolt 8% Bis-Tris Plus Gel, 17 Well	NW00087BOX
Bolt 10% Bis-Tris Plus Gel, 10 Well	NW00100BOX
Bolt 10% Bis-Tris Plus Gel, 12 Well	NW00102BOX
Bolt 10% Bis-Tris Plus Gel, 15 Well	NW00105BOX
Bolt 10% Bis-Tris Plus Gel, 17 Well	NW00107BOX
Bolt 12% Bis-Tris Plus Gel, 10 Well	NW00120BOX
Bolt 12% Bis-Tris Plus Gel, 12 Well	NW00122BOX
Bolt 12% Bis-Tris Plus Gel, 15 Well	NW00125BOX
Bolt 12% Bis-Tris Plus Gel, 17 Well	NW00127BOX
Bolt 4-12% Bis-Tris Plus Gel, 10 Well	NW04120BOX
Bolt 4-12% Bis-Tris Plus Gel, 12 Well	NW04122BOX
Bolt 4-12% Bis-Tris Plus Gel, 15 Well	NW04125BOX
Bolt 4-12% Bis-Tris Plus Gel, 17 Well	NW04127BOX

Qualifying NuPAGE Tris-Acetate and Bis-Tris mini gels

Description	Cat. no.	Description	Cat. no.
NuPAGE 7% Tris-Acetate Protein Gels, 1.0 mm, 12 Well	EA03552BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 12 Well	NP0322BOX
NuPAGE 7% Tris-Acetate Protein Gels, 1.0 mm, 15 Well	EA03555BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 15 Well	NP0323BOX
NuPAGE 7% Tris-Acetate Protein Gels, 1.0 mm, 10 Well	EA0355BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 1 Well	NP0324BOX
NuPAGE 7% Tris-Acetate Protein Gels, 1.5 mm, 15 Well	EA03585BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 2D Well	NP0326BOX
NuPAGE 7% Tris-Acetate Protein Gels, 1.5 mm, 10 Well	EA0358BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 9 Well	NP0327BOX
NuPAGE 3-8% Tris-Acetate Protein Gels, 1.0 mm, 12 Well	EA03752BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 17 Well	NP0329BOX
NuPAGE 3-8% Tris-Acetate Protein Gels, 1.0 mm, 15 Well	EA03755BOX	NuPAGE™ 4-12% Bis-Tris ZOOM™ Protein Gels, 1.0 mm, IPG Well	NP0330BOX
NuPAGE 3-8% Tris-Acetate Protein Gels, 1.0 mm, 10 Well	EA0375BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.5 mm, 10 Well	NP0335BOX
NuPAGE 3-8% Tris-Acetate Protein Gels, 1.0 mm, 2D Well	EA0376BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.5 mm, 15 Well	NP0336BOX
NuPAGE 3-8% Tris-Acetate Protein Gels, 1.5 mm, 15 Well	EA03785BOX	NuPAGE 4-12% Bis-Tris Protein Gels, 1.5 mm, 2D Well	NP0337BOX
NuPAGE 3-8% Tris-Acetate Protein Gels, 1.5 mm, 10 Well	EA0378BOX	NuPAGE 12% Bis-Tris Protein Gels, 1.0 mm, 10 Well	NP0341BOX
NuPAGE 10% Bis-Tris Protein Gels, 1.0 mm, 10 Well	NP0301BOX	NuPAGE 12% Bis-Tris Protein Gels, 1.0 mm, 12 Well	NP0342BOX
NuPAGE 10% Bis-Tris Protein Gels, 1.0 mm, 12 Well	NP0302BOX	NuPAGE 12% Bis-Tris Protein Gels, 1.0 mm, 15 Well	NP0343BOX
NuPAGE 10% Bis-Tris Protein Gels, 1.0 mm, 15 Well	NP0303BOX	NuPAGE 12% Bis-Tris Protein Gels, 1.0 mm, 1 Well	NP0344BOX
NuPAGE 10% Bis-Tris Protein Gels, 1.0 mm, 1 Well	NP0304BOX	NuPAGE 12% Bis-Tris Protein Gels, 1.0 mm, 2D Well	NP0346BOX
NuPAGE 10% Bis-Tris Protein Gels, 1.0 mm, 2D Well	NP0306BOX	NuPAGE 12% Bis-Tris Protein Gels, 1.0 mm, 17 well	NP0349BOX
NuPAGE 10% Bis-Tris Protein Gels, 1.0 mm, 9 Well	NP0307BOX		
NuPAGE 10% Bis-Tris Protein Gels, 1.5 mm, 10 Well	NP0315BOX		
NuPAGE 10% Bis-Tris Protein Gels, 1.5 mm, 15 Well	NP0316BOX		
NuPAGE 4-12% Bis-Tris Protein Gels, 1.0 mm, 10 Well	NP0321BOX		

Product to offer

Description	Cat. no.
Mini Gel Tank	A25977

* This promotion is available to customers in Albania, Bosnia & Herzegovina, Bulgaria, Croatia, Cyprus, EIRE, Estonia, Greece, Iceland, Israel, Latvia, Lithuania, Macedonia, Malta, Montenegro, Romania, Serbia, Slovenia and Turkey, who purchase five (5) boxes of qualifying Bolt mini gels or Novex WedgeWell mini gels or qualifying NuPAGE mini gels in a single purchase order at list price and include promotion code R5GEL3 for the A25977 Mini Gel Tank. Offer valid for qualifying orders received no later than June 30, 2019. Limited to one redemption per order. Offer applies to list price in effect at the time order is received. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation, or agency/institutional policy. Other restrictions may apply.

Pouring your own gels?

You can have a tank for free



Simply quote promo code **RSURE4** when ordering one of the SureCast Gel Handcast System starter packs and get a Mini Gel Tank for free†

SureCast Gel Handcast System

Pour your own polyacrylamide mini gels easily and confidently with the Invitrogen™ SureCast™ Gel Handcast System. The system features durable glass plates that are securely sealed using a single-motion load-and-lock mechanism. The unique tilt feature of the Invitrogen™ SureCast™ Gel Handcast Station minimises spillage during handcasting.

- Leak-free design—minimise the number of failed gels, waste less time
- Superior glass plate durability—up to 20 times more durable compared to other suppliers
- Unique tilt feature—minimise spillage and confidently load acrylamide solutions
- Simple component assembly—using a single-motion load-and-lock mechanism

Description	Cat. no.
SureCast™ Gel Handcast Bundle A - Hardware and Reagents	HC1000SR
SureCast™ Gel Handcast Bundle B - Hardware Only	HC1000S

Product to offer

Description	Cat. no.
Mini Gel Tank	A25977

For more information, go to thermofisher.com/proteingels

To place your order, please contact your authorised distributor

Bio Sciences Ltd - info@biosciences.ie - www.biosciences.ie



BioSciences
Delivering Excellence

† This promotion is available to customers in Albania, Bosnia & Herzegovina, Bulgaria, Croatia, Cyprus, EIRE, Estonia, Greece, Iceland, Israel, Latvia, Lithuania, Macedonia, Malta, Montenegro, Romania, Serbia, Slovenia and Turkey, who purchase one of the SureCast Gel Handcast System start packs and include promotion code RSURE4 for the A25977 Mini Gel Tank. Offer valid for qualifying orders received no later than June 30, 2019. Limited to one redemption per order. Offer applies to list price in effect at the time order is received. Cannot be combined with other discounts or promotions. Offer void where prohibited, licensed, or restricted by federal, state, provincial, or local laws or regulation, or agency/institutional policy. Other restrictions may apply.

ThermoFisher
SCIENTIFIC

For Research Use Only. Not for use in diagnostic procedures. © 2019 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **CO128532 0217**